


Phospho-Nf-Kb P65/Rela-S536 Rabbit Monoclonal Antibody For Wb, Elisa - Q04206

Item Number: CM0001649



Introduction

Phospho-specific rabbit monoclonal antibody targeting human NF-kB p65/RelA phosphorylated at Ser536. Validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for studying NF-kB signaling in inflammation, immunity, and cancer.

[Learn More](#)

Parameter	Value
Product Name	Phospho-NF-kB p65/RelA-S536 Rabbit mAb
Remarks/Alias	p65; CMCU; NFKB3; AIF3BL3; Phospho-NF-kB p65/RelA-S536
Species	Human
GeneID (Human)	5970
GeneID	5970
Immunogen	A synthetic phosphorylated peptide around S536 of human NF-kB p65/RelA (NP_068810.3)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q04206
Protein Weight	60 kDa
Shipping	Ice bag