


Phospho-Nf-Kb P65/RelA-S536 Rabbit Pab - Q04206

Item Number: CM0001489



Introduction

Phospho-specific rabbit polyclonal antibody detecting NF-kB p65 phosphorylated at serine 536. Validated for WB, IHC-P, and ELISA in human, mouse, and rat samples. Ideal for NF-kB signaling and inflammation research.

[Learn More](#)

Parameter	Details
Product Name	Phospho-NF-kB p65/RelA-S536 Rabbit pAb
Remarks/Alias	p65; CMCU; NFKB3; AIF3BL3; Phospho-NF-kB p65/RelA-S536
Species	Human
GeneID	5970 (Human)
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q04206
Protein Weight	60kDa
Shipping	Ice bag