

Phospho-Nf-Kb P65/Rela (Ser468) Rabbit Polyclonal Antibody For Wb, Elisa - Q04206

Item Number: CM0001905



Introduction

Phospho-specific rabbit polyclonal antibody targeting NF-κB p65/RelA at Ser468. Validated for WB and ELISA in human and rat samples. Ideal for studying NF-κB signaling, inflammation, and transcription regulation pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-NF-κB p65/RelA-S468 Rabbit pAb
Remarks/Alias	p65; CMCU; NFKB3; AIF3BL3; Phospho-NF-κB p65/RelA-S468
Species	Human
GeneID (Human)	5970
GeneID	5970
Immunogen	A synthetic phosphorylated peptide around S468 of human Phospho-NF-κB p65/RelA-S468 (NP_001138610.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q04206
Protein Weight	60 kDa
Shipping	Ice bag