


Phospho-Nf-Kb P65/Rela-T505 Rabbit Pab - Q04206

Item Number: CM0001625



Introduction

Phospho-specific rabbit polyclonal antibody targeting NF-kB p65/RelA phosphorylated at Thr505. Validated for western blot and ELISA using human samples. Suitable for NF-kB signaling, inflammation, and immunity research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-NF-kB p65/RelA-T505 Rabbit pAb
Remarks/Alias	p65; CMCU; NFKB3; AIF3BL3; Phospho-NF-kB p65/RelA-T505
Species	Human
GeneID (human)	5970
GeneID	5970
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q04206
Protein Weight	60kDa
Shipping	Ice bag