


Phospho-Nf-Kb2-S866 Rabbit Monoclonal Antibody For Wb, Elisa - Q00653

Item Number: CM0000964



Introduction

Phospho-specific rabbit monoclonal antibody against NF-κB2 (p100/p52) phosphorylated at S866, validated for WB and ELISA. Detects endogenous p100/p52; cross-reacts with mouse. Essential for non-canonical NF-κB pathway and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-NF-κB2-S866 Rabbit mAb
Remarks/Alias	p52; p100; H2TF1; LYT10; CVID10; LYT-10; NF-κB2; p49/p100; Phospho-NF-κB2-S866
Species	Human
GeneID (human)	4791
GeneID	4791
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around S866 of human Phospho-NF-κB2-S866 (NP_001070962.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q00653
Protein Weight	97kDa
Shipping	Ice bag