

Phospho-Nf2-S518 Rabbit Pab - P35240

Item Number: CM0000509



Introduction

Phospho-NF2-S518 Rabbit pAb detects endogenous levels of NF2/Merlin phosphorylated at serine 518. Validated for WB, IHC-P, and ELISA in human, mouse, and rat samples, supporting tumor suppressor and Hippo pathway research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-NF2-S518 Rabbit pAb
Remarks/Alias	ACN; SCH; BANF; merlin-1; Phospho-NF2-S518
Species	Human
GeneID (Human)	4771
GeneID	4771
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S518 of human NF2 (NP_000259.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35240
Protein Weight	70kDa
Shipping	Ice bag