

## Phospho-Nfkb1 (S337) Rabbit Polyclonal Antibody For Wb, Elisa - P19838

Item Number: CM0001123



### Introduction

Phospho-specific rabbit polyclonal antibody detecting NFKB1 phosphorylated at S337. Validated for Western blot and ELISA, cross-reacts with human, mouse, and rat. Ideal for NF-kappa-B signaling studies.

[Learn More](#)

| Parameter        | Value                                                                                                                       |
|------------------|-----------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Phospho-NFKB1-S337 Rabbit pAb                                                                                               |
| Remarks/Alias    | KBF1; EBP-1; NF-kB; CVID12; NF-kB1; NFKB-p50; NFKappaB; NF-kappaB; NFKB-p105; NF-kappa-B1; NF-kappabeta; Phospho-NFKB1-S337 |
| Species          | Human                                                                                                                       |
| GeneID (Human)   | 4790                                                                                                                        |
| GeneID           | 4790                                                                                                                        |
| Immunogen        | Synthetic peptide A synthetic phosphorylated peptide around S337 of human NFKB1 (NP_001158884.1).                           |
| Source           | Rabbit                                                                                                                      |
| Category         | Phosphorylated Antibodies                                                                                                   |
| Application      | WB, ELISA                                                                                                                   |
| Cross Reactivity | Human, Mouse, Rat                                                                                                           |
| SWISS            | P19838                                                                                                                      |
| Protein Weight   | 50 kDa/105 kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                                     |