


Phospho-Nfkb1-S337 Rabbit Polymab® For Wb - P19838

Item Number: CM0001039



Introduction

Rabbit monoclonal antibody against phosphorylated NFKB1 (Ser337), validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Ideal for studying NF-κB signaling in inflammation and immunity.

[Learn More](#)

Parameter	Value
Product Name	Phospho-NFKB1-S337 Rabbit PolymAb®
Remarks/Alias	KBF1; EBP-1; NF-κB; CVID12; NF-κB1; NFKB-p50; NFκappaB; NF-kappaB; NFKB-p105; NF-kappa-B1; NF-kappabeta
Species	Human
Gene ID (Human)	4790
Gene ID	4790
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S337 of human NFKB1 (NP_001158884.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P19838
Protein Weight	50kDa/105kDa
Shipping	Ice bag