

Phospho-Nfkb1-S907 Rabbit Polyclonal Antibody For Western Blot - P19838

Item Number: CM0001801



Introduction

High-quality phospho-specific NFKB1 (pS907) rabbit polyclonal antibody validated for Western blot detection of human phosphorylated NF-kappa-B p105/p50. Ideal for inflammation, immunity, and cancer signaling research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-NFKB1-S907 Rabbit pAb
Remarks/Alias	KBF1, EBP-1, NF-kB, CVID12, NF-kB1, NFKB-p50, NFkappaB, NF-kappaB, NFKB-p105, NF-kappa-B1, NF-kappabeta, Phospho-NFKB1-S907
Species	Human
Gene ID	4790
Immunogen	A phospho specific peptide corresponding to residues surrounding S907 of human NFKB1
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB
Cross Reactivity	Human
SWISS	P19838
Protein Weight	105kDa
Shipping	Ice bag