

Phospho-Nfkb1-S927 Rabbit Pab For Wb, Elisa - P19838

Item Number: CM0001388



Introduction

Rabbit polyclonal antibody specific for phosphorylated NFKB1 (Ser927). Suitable for Western blot and ELISA with human samples. A valuable tool for studying NF-kappa-B signaling in inflammation, immunity, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-NFKB1-S927 Rabbit pAb
Remarks/Alias	KBF1; EBP-1; NF-kB; CVID12; NF-kB1; NFKB-p50; Nfkapab; NF-kappaB; NFKB-p105; NF-kappa-B1; NF-kappabeta; Phospho-NFKB1-S927
Species	Human
GeneID	4790
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around S927 of human NFKB1 (NP_001158884).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P19838
Protein Weight	105kDa
Shipping	Ice bag