

Phospho-Nmdar1-S896 Rabbit Polyclonal Antibody For Wb, Elisa - Q05586

Item Number: CM0000754



Introduction

Phospho-NMDAR1-S896 Rabbit pAb specifically detects glutamate receptor ionotropic NMDA 1 phosphorylated at serine 896 in human, mouse, and rat samples. Ideal for Western blot (WB) and ELISA applications in neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-NMDAR1-S896 Rabbit pAb
Remarks/Alias	NR1; MRD8; GluN1; NMDA1; DEE101; NDHMSD; NDHMSR; NMD-R1; NMDAR1; Phospho-NMDAR1-S896
Species	Human
Gene ID (Human)	2902
Gene ID	2902
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around S896 of human NMDAR1 (NP_015566.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q05586
Protein Weight	105kDa
Shipping	Ice bag