

Phospho-P38 Mapk-T180 Rabbit Pab For Wb, Ihc-P, Elisa - Q16539

Item Number: CM0001483



Introduction

Rabbit polyclonal antibody specific for phosphorylated p38 MAPK (Thr180). Suitable for WB, IHC-P, ELISA; cross-reacts with human, mouse, rat MAPK14. Ideal for studying stress-activated kinase signaling and inflammatory pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-p38 MAPK-T180 Rabbit pAb
Remarks/Alias	RK; p38; CSBP; EXIP; Mxi2; CSBP1; CSBP2; CSPB1; PRKM14; PRKM15; SAPK2A; p38ALPHA; Phospho-p38 MAPK-T180
Species	Human
GeneID (Human)	1432
GeneID	1432
Immunogen	A synthetic phosphorylated peptide around T180 of human MAPK (NP_620581.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16539
Protein Weight	41kDa
Shipping	Ice bag