


Phospho-P38 Mapk-T180/Y182 Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q15759/P53778/O15264/Q16539

Item Number: CM0001968



Introduction

Rabbit polyclonal antibody detecting phosphorylated T180/Y182 of human p38 MAPK. Validated for WB, IF-P, IHC-P, ELISA in human, mouse, rat. Ideal for MAPK pathway and stress response research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-p38 MAPK-T180/Y182 Rabbit pAb
Remarks/Alias	MAPK14; CSBP; CSBP1; CSBP2; CSPB1; EXIP; Mxi2; PRKM14; PRKM15; RK; SAPK2A; p38; p38ALPHA; p38 MAPK; Phospho-p38 MAPK-T180/Y182
Species	Human
Gene ID (Human)	1432
Gene ID	5600/6300/5603/1432
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15759/P53778/O15264/Q16539
Protein Weight	29kDa/34kDa/35kDa/41kDa
Shipping	Ice bag