


Phospho-P38 Mapk (T180/Y182) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q16539

Item Number: CM0001283



Introduction

Phospho-p38 MAPK (T180/Y182) rabbit polyclonal antibody validated for WB, IHC-P, ELISA. Detects human, mouse, rat p38 alpha. Suitable for inflammation and stress signaling research. Q16539.

[Learn More](#)

Parameter	Value
Product Name	Phospho-p38 MAPK-T180/Y182 Rabbit pAb
Remarks/Alias	RK; p38; CSBP; EXIP; Mxi2; CSBP1; CSBP2; CSPB1; PRKM14; PRKM15; SAPK2A; p38ALPHA; Phospho-p38 MAPK-T180/Y182
Species	Human
Gene ID (Human)	1432
Immunogen	A synthetic phosphorylated peptide around T180 & Y182 of human p38 MAPK (NP_620581.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16539
Protein Weight	41kDa
Shipping	Ice bag