

Phospho-P53-S367 Rabbit Polyclonal Antibody For Wb And Elisa - P04637

Item Number: CM0001420



Introduction

Phospho-specific rabbit polyclonal antibody against p53 phosphorylated at Ser367, validated for Western blot (WB) and ELISA. Ideal for studying human p53 signaling, apoptosis, and tumor suppression in cancer research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-p53-S367 Rabbit pAb
Remarks/Alias	P53; BCC7; LFS1; BMFS5; TRP53; Phospho-p53-S367
Species	Human
GeneID (human)	7157
GeneID	7157
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around S367 of human TP53 (NP_000537.3).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P04637
Protein Weight	44kDa
Shipping	Ice bag