

Phospho-P53-S37 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P04637

Item Number: CM0001327



Introduction

Phospho-specific rabbit polyclonal antibody targeting p53 phosphorylated at Ser37. Validated for WB, IHC-P, ELISA with cross-reactivity to human, mouse, rat. Ideal for studying p53 signaling, cell cycle, and apoptosis.

[Learn More](#)

Parameter	Value
Product Name	Phospho-p53-S37 Rabbit pAb
Remarks/Alias	P53; BCC7; LFS1; BMFS5; TRP53; Phospho-p53-S37
Species	Human
Gene ID (Human)	7157
Gene ID	7157
Immunogen	A synthetic phosphorylated peptide around S37 of human p53 (NP_000537.3)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04637
Protein Weight	44kDa
Shipping	Ice bag