

Phospho-P53-S376 Rabbit Monoclonal Antibody For Wb, Elisa - P04637

Item Number: CM0001704



Introduction

Phospho-specific rabbit monoclonal antibody detecting p53 phosphorylated at serine 376. Suitable for Western blot and ELISA, reacting with human and rat p53. Ideal for apoptosis, cell cycle, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-p53-S376 Rabbit mAb
Remarks/Alias	P53; BCC7; LFS1; BMF55; TRP53; Phospho-p53-S376
Species	Human
GeneID (Human)	7157
GeneID	7157
Immunogen	Synthetic Peptide A synthetic phosphorylated peptide around S376 of human p53 (P04637).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P04637
Protein Weight	44kDa
Shipping	Ice bag