


Phospho-P53-S46 Rabbit Polyclonal Antibody For Wb, Ip, Elisa - P04637

Item Number: CM0001303



Introduction

Phospho-specific p53 (Ser46) rabbit polyclonal antibody validated for WB, IP, ELISA with human samples. Detects tumor suppressor p53 phosphorylation at Ser46, a key event in DNA damage-induced apoptosis.

[Learn More](#)

Parameter	Value
Product Name	Phospho-p53-S46 Rabbit pAb
Remarks/Alias	P53; BCC7; LFS1; BMFS5; TRP53; Phospho-p53-S46
Species	Human
Gene ID	7157
Immunogen	A synthetic phosphorylated peptide around S46 of human p53 (NP_000537.3)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human
SWISS	P04637
Protein Weight	44kDa
Shipping	Ice bag