


Phospho-P53-S6 Rabbit Pab - P04637

Item Number: CM0001479



Introduction

High-specificity Phospho-p53-S6 Rabbit Polyclonal Antibody validated for WB and ELISA, targeting human p53 phosphorylated at serine 6. Suitable for apoptosis and cell cycle studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-p53-S6 Rabbit pAb
Remarks/Alias	P53; BCC7; LFS1; BMFS5; TRP53; Phospho-p53-S6
Species	Human
GeneID (Human)	7157
GeneID	7157
Immunogen	Synthetic Peptide A synthetic phosphorylated peptide around S6 of human p53 (NP_000537.3).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P04637
Protein Weight	44kDa
Shipping	Ice bag