


Phospho-P53-S9 Rabbit Mab For Wb, Elisa - P04637

Item Number: CM0001146



Introduction

Phospho-p53-S9 Rabbit Monoclonal Antibody detects human p53 phosphorylated at serine 9, with mouse cross-reactivity. Ideal for WB and ELISA in apoptosis, cell cycle, and tumor suppressor research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-p53-S9 Rabbit mAb
Remarks/Alias	P53; BCC7; LFS1; BMF55; TRP53; Phospho-p53-S9
Species	Human
GeneID	7157
Immunogen	Synthetic Peptide - A synthetic phosphorylated peptide around S9 of human p53 (P04637)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P04637
Protein Weight	44kDa
Shipping	Ice bag