


Phospho-P53-T18 Rabbit Polyclonal Antibody For Wb, Elisa - P04637

Item Number: CM0001111



Introduction

Phospho-p53 (Thr18) rabbit polyclonal antibody validated for Western blot and ELISA. Detects human p53, a critical tumor suppressor protein involved in cell cycle arrest, apoptosis, and DNA repair.

[Learn More](#)

Parameter	Value
Product Name	Phospho-p53-T18 Rabbit pAb
Remarks/Alias	P53; BCC7; LFS1; BMFS5; TRP53; Phospho-p53-T18
Species	Human
Gene ID	7157
Immunogen	A synthetic phosphorylated peptide around T18 of human p53 (NP_000537.3)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P04637
Protein Weight	44kDa
Shipping	Ice bag