


Phospho-P70 S6 Kinase 1-T389 Rabbit Pab - P23443

Item Number: CM0001849



Introduction

Rabbit polyclonal antibody detecting endogenous phosphorylated p70 S6 Kinase 1 at T389. Validated for WB, IHC-P, ELISA with cross-reactivity to human, mouse, and rat. Ideal for mTOR pathway and cancer research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-p70 S6 Kinase 1-T389 Rabbit pAb
Remarks/Alias	S6K; PS6K; S6K1; STK14A; p70-S6K; p70 S6KA; p70-alpha; S6K-beta-1; p70(S6K)-alpha; Phospho-p70 S6 Kinase 1-T389
Species	Human
GeneID (human)	6198
GeneID	6198
Immunogen	Synthetic Peptide: A synthetic phosphorylated peptide around T389 of human P70S6K1 (NP_001258989.1?).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23443
Protein Weight	59kDa
Shipping	Ice bag