


Phospho-P70 S6 Kinase 1-T389 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P23443

Item Number: CM0000951



Introduction

Phospho-specific rabbit polyclonal antibody targeting p70 S6 Kinase 1 phosphorylated at T389, suitable for WB, IF/ICC, and ELISA applications in human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	Phospho-p70 S6 Kinase 1-T389 Rabbit pAb
Remarks/Alias	S6K; PS6K; S6K1; STK14A; p70-S6K; p70 S6KA; p70-alpha; S6K-beta-1; p70(S6K)-alpha; Phospho-p70 S6 Kinase 1-T389
Species	Human
GeneID (Human)	6198
GeneID	6198
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around T389 of human P70 S6K (NP_003152.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23443
Protein Weight	59kDa
Shipping	Ice bag