

Phospho-P70 S6 Kinase 1-T421/S424 Rabbit Mab For Wb, If/Icc, Elisa - P23443

Item Number: CM0001517



Introduction

High-quality Phospho-p70 S6 Kinase 1-T421/S424 Rabbit monoclonal antibody, validated for WB, IF/ICC, ELISA, with reactivity against human, mouse, and rat. Ideal for mTOR signaling, cell growth, and translation research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-p70 S6 Kinase 1-T421/S424 Rabbit mAb
Remarks/Alias	S6K; PS6K; S6K1; STK14A; p70-S6K; p70 S6KA; p70-alpha; S6K-beta-1; p70(S6K)-alpha; Phospho-p70 S6 Kinase 1-T421/S424
Species	Human
Gene ID (Human)	6198
Gene ID	6198
Immunogen	Synthetic Peptide A synthetic phosphorylated peptide around T421&S424 of human Phospho S6K1(NP_001258989.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23443
Protein Weight	59kDa
Shipping	Ice bag