

Phospho-P70 S6 Kinase 1-T421/S424 Rabbit Pab For Wb, Elisa - P23443

Item Number: CM0000523



Introduction

Rabbit polyclonal antibody against phosphorylated p70 S6 Kinase 1 (T421/S424). Validated for Western blot and ELISA. Suitable for human, mouse, and rat samples. Ideal for mTOR signaling and cell growth research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-p70 S6 Kinase 1-T421/S424 Rabbit pAb
Remarks/Alias	S6K; PS6K; S6K1; STK14A; p70-S6K; p70 S6KA; p70-alpha; S6K-beta-1; p70(S6K)-alpha; Phospho-p70 S6 Kinase 1-T421/S424
Species	Human
Gene ID (Human)	6198
Gene ID	6198
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around T421 & S424 of human RP56KB1 (NP_001258989.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23443
Protein Weight	59kDa
Shipping	Ice bag