


Phospho-P90Rsk (T359/S363) Rabbit Monoclonal Antibody For Wb - Q15418

Item Number: CM0001808



Introduction

Phospho-specific rabbit monoclonal antibody targeting p90Rsk/RSK1 phosphorylated at T359/S363, validated for WB, IP, IF/ICC, and ELISA applications. Reacts with human, mouse, and rat. Ideal for studying ERK/MAPK signaling, cell survival, and proliferation.

[Learn More](#)

Parameter	Value
Product Name	Phospho-p90Rsk-T359/S363 Rabbit mAb
Remarks/Alias	RSK; HU-1; RSK1; p90Rsk; MAPKAPK1; MAPKAPK1A
Species	Human
GeneID (Human)	6195
GeneID	6195
Immunogen	Synthetic Peptide: A synthetic phosphorylated peptide around T359 & S363 of human p90Rsk/RSK1/RPS6KA1 (NP_002944.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IP, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15418
Protein Weight	83kDa
Shipping	Ice bag