


Phospho-P90Rsk/Rsk1/Rps6Ka1-S363 Rabbit Pab - Q15418

Item Number: CM0001423



Introduction

Phospho-specific rabbit polyclonal antibody targeting human p90Rsk/RSK1/RPS6KA1 phosphorylated at Ser363. Validated for WB, IP, ELISA; cross-reacts with rat. Detects 83 kDa protein. For kinase signaling research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-p90Rsk/RSK1/RPS6KA1-S363 Rabbit pAb
Remarks/Alias	RSK; HU-1; RSK1; p90Rsk; MAPKAPK1; MAPKAPK1A; Phospho-p90Rsk/RSK1/RPS6KA1-S363
Species	Human
GeneID (human)	6195
GeneID	6195
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around S363 of human p90Rsk/RSK1/RPS6KA1 (NP_002944.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Rat
SWISS	Q15418
Protein Weight	83kDa
Shipping	Ice bag