


Phospho-P90Rsk/Rsk1/Rps6Ka1-S380 Rabbit Polyclonal Antibody For Wb, Ip, Elisa - Q15418

Item Number: CM0000714



Introduction

Phospho-specific rabbit polyclonal antibody targeting p90Rsk/RSK1/RPS6KA1 at Ser380. Validated for WB, IP, ELISA in human, mouse, and rat samples. Ideal for studying MAPK/ERK signaling, cell growth, and apoptosis.

[Learn More](#)

Parameter	Value
Product Name	Phospho-p90Rsk/RSK1/RPS6KA1-S380 Rabbit pAb
Remarks/Alias	RSK; HU-1; RSK1; p90Rsk; MAPKAPK1; MAPKAPK1A; Phospho-p90Rsk/RSK1/RPS6KA1-S380
Species	Human
GeneID	6195
Immunogen	A synthetic phosphorylated peptide around S380 of human p90Rsk/RSK1/RPS6KA1 (NP_002944.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15418
Protein Weight	83kDa
Shipping	Ice bag