


Phospho-P90Rsk/Rsk1/Rps6Ka1-T348 Rabbit Pab - Q15418

Item Number: CM0000916



Introduction

Phospho-specific rabbit polyclonal antibody targeting p90Rsk/RSK1 phosphorylated at Thr348. Validated for Western blot and ELISA, cross-reacts with human, mouse, and rat. Ideal for studying MAPK/ERK signaling and cell proliferation.

[Learn More](#)

Parameter	Value
Product Name	Phospho-p90Rsk/RSK1/RPS6KA1-T348 Rabbit pAb
Remarks/Alias	RSK; HU-1; RSK1; p90Rsk; MAPKAPK1; MAPKAPK1A; Phospho-p90Rsk/RSK1/RPS6KA1-T348
Species	Human
Gene ID	6195
Immunogen	A synthetic phosphorylated peptide around T348 of human p90Rsk/RSK1/RPS6KA1
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15418
Protein Weight	83kDa
Shipping	Ice bag