


Phospho-P90Rsk/Rsk1/Rps6Ka1-T359/S363 Rabbit Polyclonal Antibody For Wb, Elisa - Q15418

Item Number: CM0001165



Introduction

Phospho-p90Rsk/RSK1/RPS6KA1-T359/S363 Rabbit Polyclonal Antibody for WB and ELISA. Detects human RPS6KA1 phosphorylated at Thr359 and Ser363, a downstream effector of ERK signaling involved in cell growth and survival. Ideal for signal transduction research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-p90Rsk/RSK1/RPS6KA1-T359/S363 Rabbit pAb
Remarks/Alias	RSK; HU-1; RSK1; p90Rsk; MAPKAPK1; MAPKAPK1A; Phospho-p90Rsk/RSK1/RPS6KA1-T359/S363
Species	Human
Gene ID (Human)	6195
Gene ID	6195
Immunogen	Synthetic peptide (a synthetic phosphorylated peptide around T359 & S363 of human p90Rsk/RSK1/RPS6KA1, NP_002944.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q15418
Protein Weight	83kDa
Shipping	Ice bag