

Phospho-Pak2-S20 Rabbit Polyclonal Antibody For Wb & Elisa - Q13177

Item Number: CM0000975



Introduction

Phospho-PAK2 (S20) rabbit polyclonal antibody detects phosphorylated PAK2 in human, mouse, and rat samples by WB and ELISA. Ideal for studying PAK2 signaling, apoptosis, and cytoskeleton regulation.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PAK2-S20 Rabbit pAb
Remarks/Alias	KNO2; PAK65; PAKgamma; Phospho-PAK2-S20
Species	Human
GeneID	5062
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13177
Protein Weight	58kDa
Shipping	Ice bag