

Phospho-Pak4/Pak5/Pak6-S474/S560/S602 Rabbit Mab - O96013

Item Number: CM0000848

Introduction

Phospho-PAK4/PAK5/PAK6-S474/S560/S602 rabbit monoclonal antibody for WB and ELISA applications. Detects human, mouse, and rat targets. Specific for phosphorylated PAK family kinases in signal transduction research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PAK4/PAK5/PAK6-S474/S560/S602 Rabbit mAb
Remarks/Alias	PAK4; KIAA1142; PAK5KIAA1264; PAK6; Phospho-PAK4/PAK5/PAK6-S474/S560/S602
Species	Human
Gene ID (Human)	10298/57144/56924
Gene ID	10298/57144/56924
Immunogen	Synthetic Peptide - A synthetic phosphorylated peptide around S474 of human PAK4 (O96013)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O96013(PAK4-Antigen), Q9P286(PAK5), Q9NQUS(PAK6)
Protein Weight	64, 80, 75 kDa
Shipping	Ice bag