


Phospho-Paxillin (Y118) Rabbit Polyclonal Antibody For Wb And Elisa - P49023

Item Number: CM0000723



Introduction

Phospho-specific rabbit polyclonal antibody targeting Paxillin phosphorylated at tyrosine 118. Validated for WB and ELISA in human samples. Detects a 65 kDa band. Suitable for focal adhesion and cell migration studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-Paxillin-Y118 Rabbit pAb
Remarks/Alias	Paxillin; PXN; paxillin; Phospho-Paxillin-Y118
Species	Human
GeneID	5829
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P49023
Protein Weight	65kDa
Shipping	Ice bag