

Phospho-Pdgfr Alpha-Y1018 Rabbit Pab For Wb, Elisa - P16234

Item Number: CM0000481



Introduction

Phospho-PDGFR alpha-Y1018 rabbit polyclonal antibody validated for WB and ELISA. Detects human PDGFRA phosphorylated at Y1018. Ideal for studying PDGF signaling, cell proliferation, and receptor tyrosine kinase pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PDGFR alpha-Y1018 Rabbit pAb
Remarks/Alias	CD140A; PDGFR2; PDGFR-2; Phospho-PDGFR alpha-Y1018
Species	Human
GeneID	5156
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P16234
Protein Weight	123kDa
Shipping	Ice bag