

Phospho-Pdgfr Alpha-Y849 Rabbit Pab For Wb, Elisa - P16234

Item Number: CM0000735



Introduction

Phospho-PDGFR alpha-Y849 Rabbit pAb detects human PDGFRA phosphorylated at tyrosine 849. Validated for Western blot and ELISA, this antibody supports research into PDGF signaling, cell proliferation, and development.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PDGFR alpha-Y849 Rabbit pAb
Remarks/Alias	CD140A; PDGFR2; PDGFR-2; Phospho-PDGFR alpha-Y849
Species	Human
GeneID (Human)	5156
GeneID	5156
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around Y849 of human PDGFRA (NP_006197.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P16234
Protein Weight	123kDa
Shipping	Ice bag