


Phospho-Pdgfrb-Y740 Rabbit Pab For Wb And Elisa - P09619

Item Number: CM0001596



Introduction

Phospho-PDGFRB-Y740 Rabbit Polyclonal Antibody validated for WB and ELISA, detects PDGFRB phosphorylated at Y740 in human and mouse. Ideal for studying PDGF signaling in cell proliferation, migration, and cardiovascular development.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PDGFRB-Y740 Rabbit pAb
Remarks/Alias	IMF1; KOGS; IBGC4; JTK12; PDGFR; PENTT; CD140B; PDGFR1; PDGFR-1; Phospho-PDGFRB-Y740
Species	Human
GeneID (Human)	5159
GeneID	5159
Immunogen	Synthetic peptide - A synthetic phosphorylated peptide around Y740 of human PDGFRB (NP_002600.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P09619
Protein Weight	124kDa
Shipping	Ice bag