

Phospho-Pdha1-S293 Rabbit Mab - P08559

Item Number: CM0000421



Introduction

Phospho-PDHA1-S293 rabbit monoclonal antibody validated for WB, IP, ELISA in human, mouse, rat. Detects phosphorylated Ser293 of PDHA1, a mitochondrial pyruvate dehydrogenase subunit linking glycolysis to the TCA cycle, for metabolic disease research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PDHA1-S293 Rabbit mAb
Remarks/Alias	PDHA; PDHAD; PHE1A; E1alpha; PDHCE1A; Phospho-PDHA1-S293
Species	Human
Gene ID (Human)	5160
Gene ID	5160
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08559
Protein Weight	43 kDa
Shipping	Ice bag