


Phospho-Pdk1/Pdpk1-S241 Rabbit Mab For Wb, Elisa - O15530

Item Number: CM0000865



Introduction

This phospho-specific rabbit monoclonal antibody targets PDK1 (PDPK1) phosphorylated at serine 241. Validated for WB and ELISA, it detects human 3-phosphoinositide-dependent protein kinase 1, a key AGC kinase regulator in cell signaling.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PDK1/PDPK1-S241 Rabbit mAb
Remarks/Alias	PDK1; PDPK2; PDPK2P; PRO0461; Phospho-PDK1/PDPK1-S241
Species	Human
Gene ID (Human)	5170
Gene ID	5170
Immunogen	Synthetic peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O15530
Protein Weight	63kDa
Shipping	Ice bag