


Phospho-Pdk1/Pdpk1-S241 Rabbit Pab - O15530

Item Number: CM0000918



Introduction

Phospho-specific rabbit polyclonal antibody detects PDK1 phosphorylated at Ser241 (PDPK1, 63 kDa). Validated for WB, IF/ICC, ELISA with human, mouse, rat cross-reactivity. Ideal for insulin signaling, cancer, and kinase pathway research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PDK1/PDPK1-S241 Rabbit pAb
Remarks/Alias	PDK1; PDPK2; PDPK2P; PRO0461; Phospho-PDK1/PDPK1-S241
Species	Human
Gene ID (Human)	5170
Gene ID	5170
Immunogen	Synthetic peptide - A synthetic phosphorylated peptide around S241 of human PDK1/PDPK1 (NP_002604.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15530
Protein Weight	63kDa
Shipping	Ice bag