


Phospho-Pdk1/Pdpk1-S241 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O15530

Item Number: CM0001899



Introduction

Phospho-specific rabbit polyclonal antibody targeting PDK1/PDPK1 at Ser241. Suitable for WB, IF/ICC, ELISA in human, mouse, and rat samples. Ideal for studying insulin signaling, cell proliferation, and survival pathways.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PDK1/PDPK1-S241 Rabbit pAb
Remarks/Alias	PDK1; PDPK2; PDPK2P; PRO0461; Phospho-PDK1/PDPK1-S241
Species	Human
GeneID (Human)	5170
GeneID	5170
Immunogen	Synthetic peptide: A phospho specific peptide corresponding to residues surrounding S241 of human PDK1/PDPK1
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15530
Protein Weight	63kDa
Shipping	Ice bag