

Phospho-Pea15-S104 Rabbit Pab For Wb, Ihc-P, Elisa - Q15121

Item Number: CM0001172



Introduction

Phospho-PEA15-S104 Rabbit polyclonal antibody (pAb) validated for WB, IHC-P, ELISA. Detects human PEA15 and shows cross-reactivity with mouse and rat. Target weight ~15kDa. Ideal for apoptosis, glucose transport, and MAPK signaling research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PEA15-S104 Rabbit pAb
Remarks/Alias	PED; MAT1; HMAT1; MAT1H; PEA-15; HUMMAT1H; PED-PEA15; PED/PEA15; Phospho-PEA15-S104
Species	Human
GeneID (Human)	8682
GeneID	8682
Immunogen	A synthetic phosphorylated peptide around S104 of human PEA15 (NP_003759.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15121
Protein Weight	15kDa
Shipping	Ice bag