

Phospho-Pea15-S104 Rabbit Polyclonal Antibody For Wb, Elisa - Q15121

Item Number: CM0000859



Introduction

Phospho-PEA15-S104 rabbit polyclonal antibody (15 kDa) validated for WB and ELISA in human, mouse, and rat. Specific for the phosphorylated S104 site of PEA-15, a regulator of apoptosis and glucose transport. Ideal for diabetes and signaling research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PEA15-S104 Rabbit pAb
Remarks/Alias	PED; MAT1; HMAT1; MAT1H; PEA-15; HUMMAT1H; PED-PEA15; PED/PEA15; Phospho-PEA15-S104
Species	Human
Gene ID	8682
Immunogen	A synthetic phosphorylated peptide around S104 of human PEA15 (NP_003759.1)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15121
Protein Weight	15kDa
Shipping	Ice bag