

Phospho-Pfkfb2-S483 Rabbit Pab For Wb, Elisa - O60825

Item Number: CM0001083



Introduction

High-quality rabbit polyclonal antibody specific for phosphorylated PFKFB2 at Ser483. Validated for WB and ELISA in human samples; ideal for studying cardiac fructose-2,6-bisphosphate metabolism.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PFKFB2-S483 Rabbit pAb
Remarks/Alias	PFK-2/FBPase-2; Phospho-PFKFB2-S483
Species	Human
Gene ID (Human)	5208
Gene ID	5208
Immunogen	Synthetic peptide A synthetic phosphorylated peptide around S483 of human PFKFB2 (NP_006203.2).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O60825
Protein Weight	58kDa
Shipping	Ice bag