

Phospho-Pi3Kp85A/P55γ/P85B-Y467/Y464/Y199 Rabbit Polyclonal Antibody For Wb - P27986 / O00459 / Q92569

Item Number: CM0001669



Introduction

Phospho-specific rabbit polyclonal antibody targeting PI3K regulatory subunits p85α, p55γ, p85β (phosphorylated Y467/Y464/Y199). Validated for WB in human, mouse, rat. Ideal for PI3K signaling, cancer, and metabolic research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PI3KP85α/P55γ/P85β-Y467/Y464/Y199 Rabbit pAb
Remarks/Alias	PIK3R1/PIK3R2/PIK3R3; Phospho-PI3KP85α/P55γ/P85β-Y467/Y464/Y199
Species	Human
GeneID (Human)	5295/5296/8503
GeneID	5295/5296/8503
Immunogen	Synthetic peptide - A phospho specific peptide corresponding to residues surrounding Y467 of human PIK3R1/2/3
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB
Cross Reactivity	Human, Mouse, Rat
SWISS	P27986/O00459/Q92569
Protein Weight	42-54kDa/81-84kDa
Shipping	Ice bag