


Phospho-Pka C-Alpha (Prkaca)-T197 Rabbit Pab - P17612

Item Number: CM0001951



Introduction

Phospho-PKA C-alpha (PRKACA)-T197 Rabbit Polyclonal Antibody validated for WB and ELISA. Recognizes human, mouse, and rat PRKACA phosphorylated at Threonine 197. Suitable for kinase signaling and phosphorylation studies.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PKA C-alpha (PRKACA)-T197 Rabbit pAb
Remarks/Alias	CAFD1; PKACA; PPNAD4; Phospho-PKA C-alpha (PRKACA)-T197
Species	Human
GeneID (Human)	5566
GeneID	5566
Immunogen	Synthetic peptide: A synthetic phosphorylated peptide around T197 of human PKA C-alpha (PRKACA) (NP_002721.1).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17612
Protein Weight	41kDa
Shipping	Ice bag