

Phospho-Pka Riiα (Prkar2A)-S99 Rabbit Mab For Wb And Elisa - P13861

Item Number: CM0001208



Introduction

Phospho-specific rabbit monoclonal antibody targeting PRKAR2A at Ser99. Validated for WB and ELISA in human and rat samples. Ideal for cAMP signaling and PKA research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PKA RIIα (PRKAR2A)-S99 Rabbit mAb
Remarks/Alias	PKR2; PRKAR2; Phospho-PKA RIIα (PRKAR2A)-S99
Species	Human
GeneID (human)	5576
GeneID	5576
Immunogen	Synthetic Peptide - A synthetic phosphorylated peptide around S99 of human PRKAR2A/PKR2 (P13861).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P13861
Protein Weight	46kDa
Shipping	Ice bag