

Phospho-Pkc Alpha-S657/Y658 Rabbit Polyclonal Antibody For Wb, Elisa - P17252

Item Number: CM0001851



Introduction

Phospho-PKC alpha (S657/Y658) rabbit polyclonal antibody validated for WB and ELISA. Detects phosphorylated human PKC alpha (77 kDa), a key kinase in cell proliferation, apoptosis, and cancer biology.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PKC alpha-S657/Y658 Rabbit pAb
Remarks/Alias	AAG6; PKCA; PRKACA; PKCI +/-; PKCalpha; PKC-alpha; Phospho-PKC alpha-S657/Y658
Species	Human
GeneID (Human)	5578
GeneID	5578
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P17252
Protein Weight	77kDa
Shipping	Ice bag