


Phospho-Pkc Alpha-T410+Pkc Delta-T410 Rabbit Mab For Wb, Elisa - P17252 / Q05655

Item Number: CM0000833



Introduction

Phospho-PKC alpha-T410+PKC delta-T410 Rabbit mAb detects phosphorylated T410 of human PKC alpha and PKC delta. Validated for WB and ELISA, ideal for cell signaling and oncology research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PKC alpha-T410+PKC delta-T410 Rabbit mAb
Remarks/Alias	AAG6; PKCA; PRKACA; PKCI+/-; PKCalpha; PKC-alpha; Phospho-PKC alpha-T410+PKC delta-T410
Species	Human
GeneID (Human)	5578/5580
GeneID	5578/5580
Immunogen	Synthetic Peptide A synthetic phosphorylated peptide around T410 of human PRKCZ(NP_002735.3)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P17252/Q05655
Protein Weight	76kDa/77kDa
Shipping	Ice bag