

Phospho-Pkc Alpha-T638 + Pkc Beta II-T641 Rabbit Monoclonal Antibody For Wb, Elisa - P17252/P05771

Item Number: CM0001106



Introduction

Phospho-specific rabbit monoclonal antibody detecting PKC alpha (Thr638) and PKC beta II (Thr641). Validated for WB and ELISA with human samples. For studying PKC signaling in cancer, cardiac biology, and inflammation.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PKC alpha-T638 + PKC beta II-T641 Rabbit mAb
Remarks/Alias	PKC-beta; PKCB; PRKCB1; PRKCB2; Phospho-PKC alpha-T638 + PKC beta II-T641
Species	Human
Gene ID (Human)	5578/5579
Gene ID	5578/5579
Immunogen	A synthetic phosphorylated peptide around T638 of human PRKCA (NP_002728.2)
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P17252/P05771
Protein Weight	76kDa/77kDa
Shipping	Ice bag