


Phospho-Pkc Alpha-T638 Rabbit Mab - P17252

Item Number: CM0000945



Introduction

Anti-Phospho-PKC alpha-T638 Rabbit Monoclonal Antibody for WB, ELISA. Detects phosphorylated PKC alpha at T638 in human and mouse samples, ideal for signal transduction research.

[Learn More](#)

Parameter	Value
Product Name	Phospho-PKC alpha-T638 Rabbit mAb
Remarks/Alias	AAG6; PKCA; PRKACA; PKCI+/-; PKCalpha; PKC-alpha; Phospho-PKC alpha-T638
Species	Human
GeneID (Human)	5578
GeneID	5578
Immunogen	Synthetic Peptide A synthetic phosphorylated peptide around T638 of human PKC alpha (P17252).
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P17252
Protein Weight	77kDa
Shipping	Ice bag